



Chromacity
www.chromacitylasers.com

Chromacity Spark 1040 nm femtosecond laser system used as irradiation source for heating of hybrid gold-iron oxide nanoparticles in biological media

Purpose

Hybrid gold-iron oxide nanoparticles have been shown to hold great potential for heat triggered drug delivery [1,2]. The magnetic core of the particles enables magnetic resonance imaging whilst the gold surface coating can act as a localised heat source after laser irradiation by exploitation of the surface plasmon resonance. Previously, we have reported heating properties in agar [3], inside cells [4] and in tumour bearing xenograft cadavers [4] after irradiation with a 1064 nm laser source. The purpose of this work was to test whether the Chromacity Spark 1040 nm femtosecond laser source was A) useful for laser irradiation of these nanoparticles and B) more effective than our previous laser source.

Methodology

Agar phantoms

Hybrid nanoparticles were suspended in 2 % agar phantoms at 5 μgml^{-1} and 50 μgml^{-1} . The gels were irradiated for 60 s using the Chromacity Spark 1040 nm femtosecond system and changes in heat monitored using an Optris PI640 thermal imaging camera (Optris, Germany). The subsequent 60 s following irradiation was also monitored to understand the cooling profile of the particles. Heat dissipation was deduced from the infrared images up to 12 mm from the point of irradiation. All reported thermal change is in respect to a control agar phantom containing no nanoparticles.

In vitro

Human pancreatic adenocarcinoma (BxPC-3) cells were seeded into 6-well plates (50,000 cells/well) containing quartz coverslips (Alfa Aesar, USA) and incubated at 37 °C with 5 % CO₂. Hybrid nanoparticles (5 μgml^{-1} & 50 μgml^{-1}) were incubated with the cells for 24 h. Cells were washed with phosphate buffered saline (PBS) four times and the coverslips removed from the plate. The cells were irradiated using the Chromacity Spark 1040 nm femtosecond laser source for 60 s or 4 repetitions of 60 s with a 60 s cooling period between. After irradiation the coverslips were returned to the 6 well plate, fresh media added and incubated for 24 h. After this time the cytotoxicity was measured using trypan blue exclusion as previously described. Cytotoxicity was measured in relation to control cells with no nanoparticles.

In vivo

Female Nu/Nu mice, five weeks of age (n=3), (Charles River, UK) were kept in pathogen-free conditions (weight of mice was 20–25 g). Human pancreatic cancer cell line BxPC-3 was cultured to 90 % confluence in RPMI 1640 supplemented with 10 % fetal bovine serum and 1 % penicillin streptomycin. The cells were washed twice with cold PBS and harvested with trypsin for 10 min at 37 °C. The cells were washed three times with PBS and resuspended in 50:50 media:PBS. The tumour cell suspension (3.0 × 10⁶ cells in 100 μl of 50:50 PBS:media) was injected subcutaneously (s.c.) in the right flank of each mouse. When the tumour became palpable (approximately after one week), measurements in two dimensions with Vernier callipers were carried out twice a week and the volume of the tumours calculated according to Equation 1.

$$V = \frac{4}{3} \pi \left[\frac{(D1+D2)}{4} \right]^3 \quad (1)$$




Figure 3. Infrared monitoring of hybrid nanoparticles in tumour bearing nude mouse xenograft cadavers after 60 s irradiation with Chromacity laser source showing **A)** Thermal change over 120 s and **B)** Heat dissipation from laser source at 60 s irradiation (n=3,±SD).

Benefits

Comparing these studies with the previous reports on hybrid nanoparticle heating shows that the Chromacity Spark 1040 nm femtosecond system has major advantages compared to the previously used 1064 nm source. Consistently, we observed that the Chromacity system provides greater thermal rise both in agar and in tumours. Previously, at 5 µgml⁻¹ no notable heating was observed compared to control samples. However, in this study in xenografts, those injected with 5 µgml⁻¹ experienced 27 °C temperature increases. Additionally, there was an 8-fold increase in temperature at 50 µgml⁻¹ compared with the 1064 nm irradiation study. These exciting findings show that for biomedical application more than a 10-fold reduction in nanoparticle concentration could be administered in order to experience similar heating effects. This would not only render formulations more cost effective but also reduce any question over the toxic burden of nano-carrier or long-term accumulation. This system would enable either localised tumour ablation or thermally triggered drug delivery which is highly localised to the laser irradiation point thus minimum impact on surrounding tissue would be experienced.

Reference

1. Oluwasanmi, A., Al-Shakarchi, W., Manzur, A., Aldebasi, M. H., Elsini, R. S., Albusair, M. K., Hoskins, C. (2017). Diels Alder-mediated release of gemcitabine from hybrid nanoparticles for enhanced pancreatic cancer therapy. *J Control Release*, 266, 355-364. doi:10.1016/j.jconrel.2017.09.027
2. Malekigorji, M., Alfahad, M., Kong Thoo Lin, P., Jones, S., Curtis, A., & Hoskins, C. (2017). Thermally triggered theranostics for pancreatic cancer therapy. *Nanoscale*, 9(34), 12735-12745. doi:10.1039/c7nr02751f
3. Curtis, A. D. M., Malekigorji, M., Holman, J., Skidmore, M., & Hoskins, C. (2015). Heat Dissipation of Hybrid Iron Oxide-Gold Nanoparticles in an Agar Phantom. *Journal of Nanomedicine and Nanotechnology*. doi:10.4172/2157-7439.1000335
4. Oluwasanmi, A., Malekigorji, M., Jones, S., Curtis, A., & Hoskins, C. (2016). Potential of hybrid iron oxide-gold nanoparticles as thermal triggers for pancreatic cancer therapy. *RSC ADVANCES*, 6(97), 95044-95054. doi:10.1039/c6ra20552f

Work carried out at the School of Pharmacy, Institute of Science and Technology in Medicine, Keele University, Keele, by Dr. Clare Hoskins and the Nanopharmaceutics group.

Reproduced with thanks to Dr. Clare Hoskins and the Nanopharmaceutics group at Keele University.



Chromacity

Livingstone House, Phone: +44 (0) 131 449 4308
 Research Avenue North Email: sales@chromacitylasers.com
 Riccarton, Edinburgh Web: www.chromacitylasers.com
 EH14 4AP Twitter: @chromacityltd